



V8 MAR.2020

# Light Bond™ Sealants


## With or Without FLUORIDE

### LIGHT CURE SEALANT RESIN

**INDICATIONS FOR USE:** Light Bond™ Sealants are intended for use as a light cure sealant resin for bonding under an orthodontic adhesive.

#### **INSTRUCTIONS FOR USE:**

Dispense one to two drops of Light Bond sealant resin onto a mixing pad. With a brush, apply a thin, uniform layer of sealant to the etched enamel. It is not necessary to cure the sealant at this time; however, a cured sealant can prevent bracket floatation.

 **PRECAUTIONS:** Unpolymerized resins may cause skin sensation in susceptible persons. In case of contact with the skin, wash thoroughly with soap and water.



**TEMPERATURE LIMITATION:** Store at room temperature (15-30 °C)

**TO REORDER:** Ref# LBS3F – 3cc Light Bond Sealant w/Fluoride, LBS7F – 7cc Light Bond Sealant w/Fluoride, LBS3 – Non-Fluoride, 3cc Light Bond Sealant, LBS7 – Non-Fluoride 7cc Light Bond Sealant, LBFS7F – Light Bond Filled Sealant 7cc w/Fluoride, LBFSF – Light Bond Filled Sealant 3cc w/ Fluoride

**WARRANTY:** Reliance Orthodontic Products, Inc. recognizes its responsibility to replace products if proven to be defective. Reliance Orthodontic Products, Inc. does not accept liability for any damages or loss, direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

**Rx Only:** Federal law restricts this device to sale by or on the order of a dental professional.



1540 West Thorndale Ave., Itasca, IL 60143 U.S.A.  
1-800-323-4348 • [www.RelianceOrthodontics.com](http://www.RelianceOrthodontics.com)

Made in the USA



EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

MD

Check out our new learning  
platform Reliance U!  
[www.RelianceU.tech](http://www.RelianceU.tech)

Australian Sponsor:  
Emergo Australia  
Level 20, Tower II, Darling Park  
201 Sussex St.  
Sydney, NSW 2000 Australia



Scan for more information